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REFERENCES:

- 1- Nuttall, K.L., Klee, G.G. Analyter of haemoglobin metabolism Porphyrins, Iron and Bilirubin, Tietz Fundamentals of Clinical Chemistry, 5thEd., Burtis, C.A.& Ashwood E.R.(W.B saunders eds.Philadelphia USA), (2001),584.
- 2- Schreiber, W.E., Iron, Porphyrin and bilirubin metabolism, Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak,S.C.,(Mosby Inc.eds St Louis USA),(2003),657
- 3- Scherwin, J.E.,Liver function. Clinical Chemistry: Theory, Analysis, Correlation, 4thEd., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc.eds St Louise USA),(2003), 493 and appendix.
- 4- Tietz, N.W, Clinical guide to laboratory tests. 3th Ed., (W.B.Saunders eds. Philadelphia USA),(1995), 90.
- 5- Vassault A., et al., Protocole de validation de techniques. (Documen B, stade 3), Ann.Biol.Clin.,(1986), 44,686.
- 6- Vassault A., et al., Analyses de Biologie medicale: Specifications et norms d acceptabilite a l usage de la validation des techniques., Ann.Biol.Clin., (1997), 57 ,685.
- 7- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2th Ed., AACC Press, (1997).
- 8- Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).
- 9- Berth, M. & Delanghe, J. Protein Precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 casse reports and review of literature, Acta Clin Belg., (2004), 59, 263.