Creatinine (Jaffe-kinetic)

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$$\frac{\text{abs Sample } (\Delta\text{abs}/\text{min})}{\text{abs Calibrator } (\Delta\text{abs}/\text{min})} \times \text{Conc. Std/Cal (mg/dL)} = \text{Conc. Creatinine (mg/dL)}$$

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$$\text{Urine 24hr (mg/24 hr)} = \left[\frac{\text{Urine Creatinine (mg/dL)} \times \text{Urine Volume (ml)}}{100} \right] \times 20$$

$$\text{Urine 24hr (mg/kg/24hr)} = \left[\frac{\text{Urine Volume (ml)} \times \text{Creatinine Conc (mg/dl)}}{\text{Body Weight (kg)}} \right] \times 100 \times 20$$

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$$\frac{\text{abs Sample } (\Delta\text{abs}/\text{min})}{\text{abs Calibrator } (\Delta\text{abs}/\text{min})} \times \text{Conc. Std/Cal (mg/dL)} = \text{Conc. Creatinine (mg/dL)} \times 20$$

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
- Allston, C.A, Non-protein nitrogenous compounds and renal function. Clinical chemistry: Concepts and application, Anderson, S.C, Cockayne, S. (W.B. saunders ads. Philadelphia USA), (1993), 369.
- Newman D.J., price C.P., Non protein nitrogen metabolites. Tietz fundamentals of clinical chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. saunders ads. Philadelphia USA), (2001), 414.
- Butler, A.R., The jaffe reaction. Identification of the colored species. Clin. Chim. Acta. (1975), 59, 227.
- Vasiliades, j., Reaction of alkaline picrate with Creatinine.1.kinetics and mechanism of formation of the mono Creatinine picrate acid complex, clin. Chem., (1976), 22, 1664.
- Peter E., Gus Koerbin., Methods in Clinical chemistry, Kaplan and Pesce's : Clinical Chemistry:Theory, Analysis, Correlation, (2009), 440.
- <https://www.mayocliniclabs.com>
- Tietz, N.W., Clinical guide to laboratory tests, 4th Ed. (W.B. Saunders eds. Philadelphia USA), (2006), 316.
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI document Ep06-A, 2003.
- Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Act Clin Belg., (2004), 59, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2ndEd. AACC Press, (1997).
- Young, D.S., Effects of drugs on clinical laboratory tests, 4thEd. AACC Press, (1995).

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 Temperature limitation	 REF	Catalogue number
 Manufacture address	 Expiration date	Expiration date
 Batch code	 Date of manufacture	Date of manufacture
 In vitro diagnostic medical device	 R1	Reagent 1
 Consult instruction for use	 R2	Reagent 2