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Acetamidophenol Deoxycorticosterone (+)3,4-Methylenedioxy- Prednisolone Acetophenetidin Dextromethorphan amphetamine Prednisone N-Acetylprocainamide Diazepam (+)3,4-Methylenedioxy- Procaine Acetylsalicylic acid Diclofenac methamphetamine Promazine Aminopyrine Diflunisal Methylphenidate Promethazine Amitriptyline Digoxin Methyprylon D,L-Propranolol Amobarbital Diphenhydramine Morphine-3- D-Proproxyphene Amoxicillin Doxylamine  $\beta$ -D-glucuronide D-Pseudoephedrine Ampicillin Ecgonine hydrochloride Nalidixic acid Quinidine L-Ascorbic acid Ecgonine methylester Nalorphine Quinine D,L-Amphetamine (-)- $\psi$ -Ephedrine Naloxone Ranitidine L-Amphetamine Erythromycin Naltrexone Salicylic acid Apomorphine  $\beta$ -Estradiol Naproxen Secobarbital Aspartame Estrone-3-sulfate Niacinamide

07

Serotonin (5-Hydroxytyramine) Atropine Ethyl-p-aminobenzoate Nifedipine Sulfamethazine Benzilic acid Fenoprofen Norcodein Sulindac Benzoic acid Furosemide Norethindrone Temazepam Benzoylcegonine Genticic acid D-Norpropoxyphene Tetracycline Benzphetamine Hemoglobin Noscapine Tetrahydrocortison, Bilirubin Hydralazine D,L-Octopamine 3-Acetate ( $\pm$ )-Brompheniramine Hydrochlorothiazide Oxalic acid Tetrahydrocortison Caffeine Hydrocodone Oxazepam 3 ( $\beta$ -D-glucuronide) Cannabidiol Hydrocortison Oxolinic acid Tetrahydrozoline Chloralhydrate O-Hydroxyhippuric acid Oxycodone Thebaine Chloramphenicol 3-Hydroxytyramine Oxymetazoline Thiamine Chlordiazepoxide Ibuprofen p-Hydroxy-Thioridazine Chlorothiazide Imipramine methamphetamine D, L-Thyroxine ( $\pm$ ) Chlorpheniramine Iproniazid Papaverine Tolbutamine Chlorpromazine ( $\pm$ ) - Isoproterenol Penicillin-G Triamterene Chlorquine Isoxsuprine Pentazocine Trifluoperazine Cholesterol Ketamine Pentobarbital Trimethoprim Clomipramine Ketoprofen Perphenazine Trimipramine Clonidine Labetalol Phencyclidine Tryptamine Cocaine hydrochloride Levorphanol Phenelzine D, L-Tryptophan Codeine Loperamide Phenobarbital Tyramine Cortisone Maprotiline Phentermine D, L-Tyrosine (-) Cotinine Meprobamate L-Phenylephrine Uric acid Creatinine Methadone  $\beta$ -Phenylethylamine Verapamil Methoxyphenamine Phenylpropanolamine Zomepirac

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Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986 2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

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Method	Hana Gene Result		Total Result	
	Results	Positive		Negative
The THC Rapid Test Dipstick	Positive	89	2	91
	Negative	1	108	109
Total Result		90	110	200
Sensitivity: > 98.8% (93.9% ~ 99.9%)		Specificity: 98.1% (93.5% ~ 99.7%)		
Accuracy: 99.5% (95.6% ~ 99.6%)		95% Confidence Intervals		

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Method	Rapid Test		Total Result	
	Results	Positive		Negative
The THC Rapid Test Dipstick	Positive	91	5	96
	Negative	4	100	104
Total Result		95	105	200
Sensitivity: 95.7% (89.5% ~ 98.8%)		Specificity: 95.2% (89.2% ~ 98.4%)		
Accuracy: 95.5% (91.6% ~ 97.9%)		95% Confidence Intervals		

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11-nor- $\Delta^9$ -THC-9 COOH Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
25	-50%	30	30	0
37.5	-25%	30	21	9
50	Cut-off	30	9	21
62.5	+25%	30	1	29
75	+50%	30	0	30
150	3X	30	0	30

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